

HIGH-PRECISION CONTINUOUS AND REAL-TIME MEASUREMENT OF O₂ USING CAVITY RING-DOWN SPECTROSCOPY FOR PHOTOSYNTHETIC LIGHT AND CO₂ RESPONSE CURVES

Magdalena Hofmann¹, David Kim-Hak¹, John Hoffnagle¹, Doug Lynch², Mark Johnson²

¹Picarro, Inc., 3105 Patrick Henry Drive, Santa Clara, California, USA

²LI-COR Biosciences, Inc, 4647 Superior Street, Lincoln, Nebraska, 68504, USA



1. ABSTRACT

Photosynthesis is a complex process that consumes carbon dioxide and water to produce oxygen and glucose. Studies to investigate leaf-level photosynthetic activity have been conducted using systems that control certain parameters such as light intensity or partial pressure of CO₂. For instance, the LI-COR LI-6800 Portable Photosynthesis System enables the user to control and measure light, temperature, CO₂ and water vapor concentrations. Measurement of O₂ fluxes concurrent with CO₂ exchange and other parameters is a technical challenge since a very high level of precision on a large O₂ background is required. Only few analytical methods including mass spectrometry, fuel cell, ultraviolet and paramagnetic cells can achieve the required high-precision O₂ analysis.

Here we present new developments of a high-precision gas analyzer that utilizes the technique of Cavity Ring-Down Spectroscopy to measure oxygen concentrations. Its compact and rugged design combined with high-precision and long-term stability allow the user to deploy the instrument in the field for continuous monitoring of atmospheric oxygen level. Measurements have a 1-σ 5-minute averaging precision of 1-2 ppm for O₂ over a dynamic range of 0-25%. We present collaborative work with LI-COR where we coupled the Picarro G2207-i O₂ analyzer with the LI-6800 Portable Photosynthesis System to enable O₂ analysis for laboratory studies of photosynthesis. To validate the setup, we conducted light and CO₂ response experiments.


2. INSTRUMENTATIONS

CRDS Technology and Picarro G2207-i Analyzer

Cavity Ring-Down Spectroscopy (CRDS) Technology utilizes the unique infrared absorption spectrum of gas-phase molecules to quantify the concentration of (and sometimes isotopes of) H₂O, CO₂, CH₄, N₂O, CH₂O, NH₃.

CRDS Features:

- Small 3-mirrored cavity ~ 35 cc
- Long effective path-length (> 10 km)
- Time-based measurement
- Laser is switched on and off, and scanned across wavelengths



The Picarro G2207-i analyzer uses the CRDS technology to measure continuously δ¹⁸O in O₂ and O₂ concentrations.

G2207-i Performance Specifications:

[O ₂] Mode	
Precision, dry [O ₂] at ambient concentration (1-σ, 5 sec/5 min, at 21% O ₂)	<20 ppm/<2 ppm
Max Drift at STP O ₂ (over 24 hrs, peak-to-peak, 1 hr interval average, at 21% O ₂)	<6 ppm
[O ₂] Operating Range	5–25%
Precision [H ₂ O] (1-σ, 5 sec)	5 ppm + 0.1% of reading

LI-6800 Portable Photosynthesis System



The LI-COR LI-6800 Portable Photosynthesis System probes the carbon fixation reactions of photosynthesis by calculating fluxes for CO₂ (assimilation) from the difference in measured concentrations of CO₂ and H₂O entering and exiting a leaf cuvette.

- Flow control
- CO₂ control
- Water vapor control
- Temperature control
- CO₂ analysis (IRGA)
- H₂O analysis (IRGA)
- Fluorometer



3. METHOD


Flow configurations




In order to reduce the uncertainty in CO₂ and O₂ gas exchange fluxes it is best to operate the setup with a large leaf area (6x6cm chamber) and low flow rate.

4. RESULTS

Light-response and CO₂ Response curves



Electron transport rate (ETR) for Zea mays



5. CONCLUSION

- High-precision continuous and real-time measurement of O₂ with the Picarro G2207-i (<2ppm 1-σ uncertainty with 300 seconds averaging)
- Straightforward coupling of the G2207-i O₂ analyzer and the LI-COR LI-6800 Portable Photosynthesis System for combined O₂ and CO₂ measurements for leaf-level physiological research
- Light response curves for the C₃ plant *Phaseolus vulgaris* (bean) and CO₂ response curve for the C₄ plant *Zea mays* (corn) are in line with expectations (expected stoichiometry for CO₂:O₂ is about 1, see [1, 2])
- Electron transport rates (ETR) for *Zea mays* from fluorescence measurements and derived from CO₂ and O₂ gross fluxes agree well.

Interested in learning more?

- Contact Magdalena Hofmann (mhofmann@picarro.com)
- Visit www.picarro.com.
- Visit www.licor.com/env

REFERENCES

- [1] Von Caemmerer. Biochemical Models of Leaf Photosynthesis. CSIRO publishing, 2000.
 [2] Ruuska, Badger, Andrews, and Von Caemmerer, Photosynthetic electron sinks in transgenic tobacco with reduced amounts of rubisco: Little evidence for significant Mehler reaction, *Journal of Experimental Botany* 51, 2000: 357-68