

# Method for flux and isotopic measurements from soil chambers using cavity ring-down spectroscopy (CRDS)

Gregor Lucic<sup>1</sup>, David Kim-Hak<sup>1</sup>, Renato Winkler<sup>1</sup>

<sup>1</sup>- Picarro Inc., 3105 Patrick Henry Drive, Santa Clara, California, USA

## 1. INTRODUCTION

- Cavity Ring-Down Spectroscopy (CRDS) is a powerful tool for the determination of soil fluxes because it enables rapid, concurrent, high-precision measurements of multiple compounds and their isotopologues.
- Quantifying the flux of key elements (C, O, H, N) from a variety of sources and sinks has never been easier, due to the portability and modularity of the CRDS-Soil Chamber methodology.
- We present a series of methodology and research highlights that demonstrates how these system may be deployed in the field.

## 2. MODULARITY

### i Picarro Analyzer

The G2000, G4000 and G5000 family of Picarro (CRDS) analyzers may be used to collect concentrations of multiple gas species (CO<sub>2</sub>, CH<sub>4</sub>, N<sub>2</sub>O, NH<sub>3</sub>, H<sub>2</sub>O, and others) and isotopologues (iCO<sub>2</sub>, iCH<sub>4</sub>, iH<sub>2</sub>O, iN<sub>2</sub>O and iO<sub>2</sub>). In many cases at the same time, with high levels of precisions, low drift and strong interference corrections.



### ii Commercial or Homemade Chambers

Chamber selection is not limited to commercial brands (LI-COR, Eosense), but may entail the use of custom-made chambers of variable complexities. With either, multiple-chamber deployments via multiplexers are supported.




### iii Application Specific Add-ons

To support field deployments, other modules or accessories may be included to enhance the system. GPS antenna, anemometer data, battery power, vehicles support, additional filters, desiccants, etc.


## 3. CONFIGURATIONS

### Option A - Parallel



For more, see LI-COR App Note 138

### Option B - In-line




For more, see Eosense App Note AN0003

## 4. ADV. APPLICATION SUPPORT

### Water and Particle Removal

Water traps, Acrodisc®, Nafion®, Drierite™ and other desiccants



### Discrete Sampling

Collect and analyze discrete, 20 ml headspace samples for concentration and/or isotopic data.




### Mobility

GPS integration, battery and generator support



### Advanced Gas Chemistry

Custom carrier gases and isotopic enrichment studies.



## 5. CASE STUDY – TEXAS A&M

**Objective:** Investigate the effect of agricultural management practices (tillage and cover cropping) on greenhouse gas emissions from soil in corn, sorghum and soybean organic cropping systems in Texas



In this setup, eight Li-COR long-term chambers are connected to a Picarro G2508 analyzer housed inside a trailer in the field




Figure 1 - Half-hourly soil CO<sub>2</sub> flux measured in the 2017 growing season (April – July) from corn, sorghum and soybean plots under conventional tillage with and without cover crops.

### Summary

- In 2017, significant differences in soil CO<sub>2</sub> flux were observed between treatments; however m 3 minute measurement time that we used per chamber was inadequate to capture N<sub>2</sub>O and CH<sub>4</sub> fluxes
- In 2018, measurement time/chamber was increased to 10 minutes which is sufficient to capture N<sub>2</sub>O and CH<sub>4</sub> fluxes from the treatment plots




Figure 2 - N<sub>2</sub>O flux from a single chamber.

For more information please contact Prof. Nithya Rajan, Dept. of Soil and Crop Sciences, Texas A&M University, E-mail: [nrajan@tamu.edu](mailto:nrajan@tamu.edu)

## 6. SUMMARY AND CONCLUSIONS

- The pairing of CRDS technology with commercial and/or homemade flux solutions opens doors to new types of research and datasets.
- The portability of this technology, paired with its ability to make high-precision, multiple-species and isotopic measurements at the same time will enrich present and future soil flux studies.

### CONTACT DETAILS

Gregor Lucic (Ph.D), Application Scientist, Picarro Inc.

[glucic@picarro.com](mailto:glucic@picarro.com)