Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





- 1 High-precision atmospheric oxygen measurement comparisons between a newly built
- 2 CRDS analyzer and existing measurement techniques

3

- 4 Tesfaye A. Berhanu<sup>1,2</sup>, John Hoffnagle<sup>2</sup>, Chris Rella<sup>2</sup>, David Kimhak<sup>2</sup>, Peter Nyfeler<sup>1</sup>, Markus
- 5 Leuenberger<sup>1</sup>
- <sup>6</sup> Climate and Environmental Physics, Physics Institute and Oeschger Centre for Climate Change Research,
- 7 University of Bern, Bern, Switzerland
- 8 <sup>2</sup>Picarro Inc., 3105 Patrick Henry Drive, Santa Clara, CA, USA

9

10

11

12

13

14

15

16

17

18

19

20

21

22

23

24

#### Abstract

Carbon dioxide and oxygen are tightly coupled in land-biospheres CO<sub>2</sub> - O<sub>2</sub> exchange processes, while they are not coupled in oceanic exchange. For this reason, atmospheric oxygen measurements can be used to constrain the global carbon cycle, especially oceanic uptake. However, accurately quantifying the small (~1-100 ppm) variations in O<sub>2</sub> is analytically challenging due to the very large atmospheric background which constitutes about 20.9 % (~209500 ppm) of atmospheric air. Here we present comprehensive laboratory and field studies for a newly developed high-precision oxygen mixing ratio and isotopic composition analyzer (Picarro G-2207) that is based on cavity ring-down spectroscopy (CRDS). From the laboratory tests, we have calculated a short-term precision (standard error of one-minute measurements) of < 1 ppm for this analyzer based on measurements of eight standard gases analyzed for two hours consecutively. In contrast to the currently existing techniques, the instrument has an excellent long-term stability and therefore a calibration every 12 hours is sufficient to get an overall uncertainty of < 5 ppm. Measurements of ambient air were also conducted at the High-Altitude Research Station, Jungfraujoch and the

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.



25

26

27

28

29

30

31

32

33

34

35

36

37

38

39

40

41

42

43

44

45

46

47

48



Beromünster tall tower in Switzerland. At both sites, we observed opposing and diurnally varying  $CO_2$  and  $O_2$  profiles due to different processes such as combustion, photosynthesis and respiration. Based on the combined measurements at Beromünster tower, we determined height dependent  $O_2$ : $CO_2$  oxidation ratios varying between -0.98 to -1.60, which increase with the height of the tower inlet, possibly due to different source contribution such as natural gas combustion with high oxidation ratio and biological processes which are at the lower end.

### 1. Introduction

Atmospheric oxygen comprises about 20.9 % of the global atmosphere and in the past decade its concentration decreased at a rate of ~ 20 per meg yr<sup>-1</sup> (Keeling and Manning, 2014) mainly associated with the increase in fossil fuel combustion. In contrast, the global average atmospheric CO<sub>2</sub> mixing ratio increased to 402.8 ppm averaged over 2016 (predicted to grow by 2 % in 2017) since its preindustrial value of 280 ppm (Le Quéré et al., 2017). As the variability of atmospheric oxygen is directly linked to the carbon cycle, both its short and long-term observations can be used to better constrain the carbon cycle. For example, since first suggested by Keeling and Shertz (1992) the long-term trends derived from concurrent measurements of atmospheric CO<sub>2</sub> and O<sub>2</sub> have been widely used to quantify the partitioning of atmospheric CO<sub>2</sub> between the land-biosphere and oceanic sinks (Battle et al., 2000; Goto et al., 2017; Manning and Keeling, 2006; Valentino et al., 2008). This method hinges on the linear coupling between CO2 and O2 with an oxidation ratio (OR) of 1.1 for the terrestrial biosphere photosynthesis-respiration processes ( $\alpha_b$ ) and 1.4 for fossil fuel combustion ( $\alpha_f$ ) while they are decoupled for oceanic processes. Meanwhile, the short-term variability in atmospheric oxygen can be used to estimate marine biological productivity and air-sea gas exchange (Keeling et al., 1998; Nevison et al., 2012). However, the accuracy of these estimates is primarily linked to the accuracy and precision of atmospheric O<sub>2</sub> measurements

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





and the assumed ORs for the different processes which are highly variable in contrast to atmospheric  $CO_2$  that can be well measured within the precision guidelines set by the Global Atmospheric Watch (GAW) ( $\pm$  0.1 ppm for the northern hemisphere).

Currently there are several techniques that can measure atmospheric O<sub>2</sub> variations as oxygen concentration based on interferometric, paramagnetic, UV absorption and fuel cell technology (Keeling, 1988a; Manning et al., 1999; Stephens et al., 2007) or as O<sub>2</sub>/N<sub>2</sub> ratios to account for the large background effect using Gas chromatography with thermal conductivity detector (GC-TCD) or Gas chromatography coupled to mass spectrometry (GC-MS) (Bender et al., 1994; Tohjima, 2000). Despite the fact that these techniques have been commercially available for more than two decades, accurate quantification of atmospheric oxygen variability remains challenging primarily because the small ppm-level atmospheric oxygen signal rides on a ~ 210,000 ppm background, which places stringent requirements on the precision and drift of the analysis methods. The techniques listed above struggle to routinely achieve the necessary performance for various reasons, including i) instability over time that requires frequent measurement interruption for calibration, ii) measurement bias with ambient and sample temperature and/or pressure, and/or iii) systematic errors in the measurement due to other atmospheric species. Further, some techniques require the use of consumables and rely on high vacuum, which complicates field deployment.

In this manuscript we describe a new high precision oxygen concentration and isotopic composition analyzer by Picarro Inc., Santa Clara, USA (G-2207) based on CRDS technology. Here, we will introduce the analyzer design principles in details, describe the unique features of the analyzer and evaluate its performance based on various independent laboratory and field tests by comparing it with currently existing techniques. Then, we will

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.



74

75

76

77

78

79

80

81

82

83

84

85

86

87

88

89

90

91

92

93

94

95



present and interpret our observations based on field measurements. Finally, we will conclude its overall performance and provide recommendations and possible improvements.

### 2. Analyzer design principles

The analyzer described here is derived from the Picarro G2000 series of CRDS analyzers. The basic elements have been described elsewhere (Crosson, 2008; Martin et al., 2016; Steig et al., 2014): briefly, the instrument is built around a high-finesse, traveling-wave optical cavity, which is coupled to either of two single-frequency DFB-stabilized semiconductor lasers. One cavity mirror is mounted on a piezoelectric translator (PZT) to allow fine tuning of the cavity resonance frequencies. A semiconductor optical amplifier between the laser sources and the cavity boosts the laser power and serves as a fast-optical switch. The cavity body is constructed of invar and enclosed in a temperature stabilized box  $(T = 45^{\circ} C, \text{ stabilized to approximately } 0.01 {\circ} C)$  for dimensional and spectroscopic stability. A vacuum pump pulls the gas to be sampled through the cavity and a proportional valve between the cavity and the pump maintains the sample pressure in the cavity at a value of 340 hPa, with variations on the order of 1 Pa. The instrument has a wavelength monitor, based upon measurements of interference fringes from a solid etalon, which is used to control the laser wavelength by adjusting the laser temperature and current. A high-speed photodiode monitors the optical power emerging from the cavity. The instrument's data acquisition system sweeps the laser frequency over the spectral feature to be measured, modulates the laser output to initiate ring-downs, and fits the ring-down signal to an exponential function to generate a spectrogram of optical loss versus laser frequency. Subsequent program modules compare the measured loss spectrum to a spectral model, using non-linear least-squares fitting (Press et al., 1986) to find the best-fit model parameters and thereby obtain a quantitative measure of the absorption due to the target molecule, and finally apply a calibration factor to

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





the optical absorption to deduce the molecular concentration. When operating in its normal gas analysis mode, the instrument acquires about 200-300 ring-downs per second and achieves a noise equivalent absorption of typically about 10<sup>-11</sup> cm<sup>-1</sup> Hz<sup>-1/2</sup>, with some variation between instruments.

The primary goal when designing this analyzer was to measure the molecular oxygen concentration with few-per-meg level precision and stability. In this context operational stability is as important as signal-to-noise. Our experience has been that the most stable operation of the analyzer is achieved when the optical phase length of the cavity is held as nearly constant as possible. In this case the free spectral range (FSR) of the temperature stabilized, invar ring-down cavity provides a better optical frequency standard than the etalon-based wavelength monitor, which in turn allows more consistent measurements of absorption line width and integrated absorption line intensity (Steig et al., 2014). For a small, field-deployable instrument, it is not practical to stabilize the absolute frequencies of the cavity modes to an optical frequency standard (Hodges et al., 2004) but the oxygen lines themselves, under conditions of constant temperature and pressure, provide an adequate frequency reference. The oxygen spectrum was also used to calibrate the FSR, by comparing a wide (approximately 10 cm<sup>-1</sup>) FSR-spaced spectrum with the Hitran database (Rothman et al., 2013).

To determine molecular oxygen concentration, the analyzer measures absorption of the Q13Q13 component of the  $a^1\Delta_g \leftarrow X^3\Sigma_g^-$  band, at a frequency of 7878.805547 cm<sup>-1</sup>, according to the latest edition of Hitran (Gordon et al., 2017). This is one of the strongest near-infrared lines of oxygen, well separated from other oxygen lines, and reasonably free of spectral interference from water, carbon dioxide, methane, and other constituents of clean air. The spectral model for this line was developed using reference spectra of clean, dry, synthetic

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.



120

121

122

123

124

125

126

127

128

129

130

131

132

133

134

135

136

137

138

139

140

141

142

143



air that were acquired with the same hardware as in the field-deployable analyzer, but with special-purpose software that allows it to operate as a more general spectrometer.

Recently, considerable work has been done to advance the understanding of spectral line shapes and to define functional representations that better describe the processes that determine spectral line shapes than does the Voigt model (Hartmann et al., 2008; Tennyson et al., 2014). Line shape studies have been published for the 1.27 µm band of O<sub>2</sub> (Fleisher et al., 2015; Lamouroux et al., 2014), though not to our knowledge for the Q branch. The apparatus used here is not capable of spectroscopic studies of comparable precision; the absolute temperature and pressure monitoring and especially the frequency metrology are far too crude for that purpose. Our goal is merely to define a simple model of the Q13Q13 line that is adequate for least-squares retrievals of the O<sub>2</sub> absorption under the limited range of conditions (stabilized temperature and pressure) that the operational analyzer experiences in the field. The CRDS analyzers use the Galatry function (Varghese and Hanson, 1984), which is distinctly better than the Voigt and still easily and quickly evaluated for line shape modeling. Ultimately, the usefulness of the spectral model is to be evaluated by the precision and stability of the O<sub>2</sub> measurements when compared with established techniques. We also note at this point that Sironneau, Fleisher, and Hodges have made detailed measurements of lines in the R branch of the  $a^1\Delta_g \leftarrow X^3\Sigma_g^-$  band and observed departures from simple, linear absorption, which they interpret as arising from collision-induced absorption (Fleisher et al., 2015). This has two important consequences for O<sub>2</sub> monitoring: the line strength is not independent of sample pressure, and optical absorption is not linear in laser intensity. We do not expect these effects to be too severe for our application because the ring-down cavity is stabilized to a very narrow range of temperature and pressure. In addition, the optical power in the ringdown cavity set by the ring-down detector threshold, which is used to trigger the

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





laser shutoff and subsequent ring-down waveform acquisition. The fact that all ring-downs occur at the same intracavity power should minimize the effect of collision-induced absorption. We have observed some excess noise on the ring-down time constants for the highest loss points at the peak of the Q13Q13 line, which might have to do with the fitting of the ring-down signal if absorption is not linear, but we cannot be certain of this explanation at present.

For spectral model development, this spectrometer has the drawback that the cavity FSR, equal to about 0.0206 cm<sup>-1</sup>, is too large to reveal much detail of the absorption line shape, even with the simplifying assumption of a Galatry line shape. We therefore acquired a set of four interleaved spectra, with the PZT-actuated mirror moved to offset the cavity modes of the individual FSR-spaced spectra by one-fourth of an FSR. The precise offsets were determined from fits to the strong and well-isolated O<sub>2</sub> lines in the spectra. From the consistency of the fitted line centers, we estimate that the positioning of the interleaved spectra was accurate to approximately 10 MHz. The spectrum of the Q13Q13 line acquired in this manner is shown in Figure 1, together with the best-fit Galatry function. It stands out that the residuals that are largely an odd function of detuning from the line center: this shows the limitations of the Galatry model in this case, since the Galatry function is purely even about the line center. The shape of the absorption line in this model is specified by two dimensionless parameters: the collisional broadening parameter

163 
$$y = \gamma / \sigma_D$$
 (1)

and the collisional narrowing parameter

$$z = \beta / \sigma_D \tag{2}$$

where γ is the frequency of broadening transitions, β is the frequency of narrowing collisions,
and σ<sub>D</sub> is the Doppler width of the transition, given by

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





 $\sigma_D = v_0 (2k_B T/Mc^2)^{1/2}$ 168 (3) 169 where v<sub>0</sub> is the transition frequency, k<sub>B</sub> is Boltzmann's constant, T is the sample temperature, 170 M is the molecular mass, and c is the speed of light. Figure 2 shows the values of y and z obtained from spectra acquired in the same way as Figure 1, as a function of cavity pressure. 171 The values depend linearly on pressure, as expected from the Galatry model, but the 172 unconstrained linear fits do not go precisely through the origin. It is not clear whether this 173 represents a breakdown of the Galatry model or simply reflects the limited quality of the data 174 set. The slope of y can be converted to an air-broadened collisional width  $\gamma_{air} = 0.0442 \text{ cm}^{-1}$ 175 <sup>1</sup>/atm, which agrees with the Hitran value of 0.0460 cm<sup>-1</sup>/atm (Rothman et al., 2013) to within 176 177 the uncertainty estimate stated by Hitran. The slope of z can be interpreted in terms of the optical diffusion coefficient (Fleisher et al., 2015), yielding D = 0.285 cm<sup>2</sup> s<sup>-1</sup>, compared to 178 the literature value of 0.233 cm<sup>2</sup> s<sup>-1</sup> for O<sub>2</sub> in air at 45 °C (Marrero and Mason, 1972). 179 Although the anticipated use of the analyzer is for ambient air samples having a very small 180 range of O<sub>2</sub> concentrations, we did investigate the variation of the line shape in binary 181 mixtures of O<sub>2</sub> and N<sub>2</sub> shown in Figure 3. The error bars are taken from the output of the 182 Levenberg-Marquardt fitting routine (Press et al., 1992). The dependence of the collisional 183 broadening parameter z on O<sub>2</sub> mole fraction was considered too small to be significant, but 184 the variation in y was used in the subsequent analysis of the air samples. Note that Wójtewicz 185 et al. (Wójtewicz et al., 2014) also found collisional broadening coefficients for nitrogen to be 186 slightly larger than for oxygen in measurements of one O<sub>2</sub> line in the B-band. 187 188 The primary goal in designing the analyzer was to achieve high enough precision to 189 make meaningful measurements of O2 in clean atmospheric samples. Although the current 190 best practice for such high-precision measurements is to work with dried samples, we decided

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.



191

192

193

194

195

196

197

198

199

200

201

202

203

204

205

206

207

208

209

210

211

212

213

214



to include high precision measurements of water vapor. There were two reasons for this decision: one is to serve as a monitor for residual water vapor, which is difficult to remove completely from the ring-down cavity and associated sample handling hardware, and the second and more ambitious reason was to see how well the effect of water vapor could be corrected for in measurements of undried ambient air. While it was considered unlikely that measurements of undried air could compete in accuracy with those of dried air, it might be possible to correct for water vapor well enough to enable useful measurements in some circumstances without the expense and inconvenience of drying the sample. For this purpose, a second laser was added, which probes the  $7_{1,6} \rightarrow 8_{4,5}$  component of the 2v<sub>3</sub> band of water vapor, at a frequency of 7816.75210 cm<sup>-1</sup> (Gordon et al., 2017). The Galatry model was used to fit spectra of synthetic air humidified to various levels of water vapor concentration. These fits also included two other nearby, very weak water lines, with intensities less than 1% of the strong transition, in order that their absorption should not perturb the line shape of the main transition. Results for the shape of the 7816.75210 cm<sup>-1</sup> line are shown in Figure 4. At the level that we can measure, only the y-parameter has a meaningful variation with water concentration. From the linear fit one obtains a pressure broadening coefficient for air,  $\gamma_{air} =$  $0.0752 \text{ cm}^{-1}/\text{atm}$ , in reasonable agreement with the Hitran value  $\gamma_{air} = 0.0787 \text{ cm}^{-1}/\text{atm}$ (Gordon et al., 2017), and a self-broadening coefficient  $\gamma_{\text{self}} = 0.413 \text{ cm}^{-1}/\text{atm}$ , to be compared with the Hitran value  $\gamma_{self} = 0.366$  cm<sup>-1</sup>/atm. Since the uncertainty estimate for the Hitran values is 10 % to 20 %, this level of agreement seems reasonable. We also looked at absorption from water near the Q13Q13 absorption line of O2. These spectra were measured in a background of pure nitrogen to reveal the very weak lines interfering with the O<sub>2</sub> measurement. Without the strong O<sub>2</sub> lines, it was impossible to interleave FSR-spaced spectra, so in this case the frequency axis comes from the analyzer's

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





wavelength monitor. The upper panel of Figure 5 shows the spectrum of saturated water vapor in nitrogen, together with a fit to a Voigt model of the molecular lines. The measurement was made at a pressure of 340 hPa and temperature of 45° C. The two most prominent features in this spectrum are actually the Q17R16 and Q13Q13 lines from traces of O<sub>2</sub> remaining in the sample while the other features are from water. The lower panel of Figure 5 shows the lines tabulated in Hitran. Immediately after the data in Figure 5 were acquired, measurements were also made at 7816.85210 cm<sup>-1</sup>, to establish the relationship between the absorption strengths in the two spectral regions. All the water lines that were observed, in both spectral regions, are from the dominant 161 isotopolgue of water, so changes in isotopic composition of atmospheric water does not lead to variation in the relative strengths of the lines we measure. Hitran simulations for molecules other than water that are expected to be present in clean, ambient air indicate that direct interference with the Q13Q13 line should be negligible at the level of precision considered here. In the case of CO<sub>2</sub>, the dilution of oxygen due to 400 ppm of CO<sub>2</sub> is significant, and larger than any direct spectral interference.

Finally, we investigated the influence of water vapor on the shape of the O<sub>2</sub> Q13Q13 line. Switching between the two lasers sources, we acquired FSR-spaced spectra of humidified synthetic air, alternately covering the 7817 cm<sup>-1</sup> and 7878 cm<sup>-1</sup> regions. Individual spectra were acquired in less than 2 s, so changes in water vapor concentration between spectra were small. These spectra, with frequency resolution of 0.0206 cm<sup>-1</sup>, were analyzed by nonlinear least-squares fitting with the following spectral models: the 7817 cm<sup>-1</sup> spectra were modeled as the sum of an empty-cavity baseline having an adjustable offset level and slope and three water peaks and the two weak perturbing peaks. The molecular absorption of the main peak was expressed as an adjustable amplitude, A<sub>w</sub>, multiplying a dimensionless, areanormalized Galatry function (Varghese and Hanson, 1984). The weak perturbers were

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.



239

240

241

242

243

244

245

246

247

248

249

250

251

252

253

254

255

256

257

258

259

260

261

262



modeled by Voigt profiles with amplitudes and line widths that constrained to be in fixed ratios to the strong line, and therefore added no new degrees of freedom to the fitting procedure. Since the amplitude A<sub>w</sub> multiplies an area-normalized shape function, it is essentially equivalent to the area of the absorption line, to the extent that the Galatry model provides a valid description of the line shape. The Doppler width of the Galatry function was fixed based on the measured cell temperature, the y-parameter was allowed to vary, and the zparameter was constrained to be proportional to y, based on the earlier measurements. In addition, the center frequency of the Galatry function was adjusted to match the data set, giving a total of five free parameters for this fit. The 7878 cm<sup>-1</sup> spectra were modeled with an adjustable baseline offset and slope and molecular absorption amplitude, A<sub>02</sub>, describing the Q13Q13 O<sub>2</sub> line. Here, too, the y-parameter and centration of the O<sub>2</sub> lines were allowed to adjust, and the z-parameter was constrained to be proportional to y. The weak water lines interfering with oxygen absorption were included in the model, but with no additional free parameters, rather the amplitudes were preset based on the measured water absorption at 7817 cm<sup>-1</sup> and the previously determined amplitude relationships between the water lines. Collisional broadening of the Q13Q13 O<sub>2</sub> line by water vapor is shown in Figure 6. From the linear fit one obtains a coefficient for collisional broadening of the Q13Q13 line by water vapor of  $\gamma_{\text{water}} = 0.0442 \text{ cm}^{-1}/\text{atm}$ . We are not aware of previous measurements of this quantity. The alternating measurements at 7817 cm<sup>-1</sup> and 7878 cm<sup>-1</sup> also calibrated the relationship between water mole fraction and the absorption at 7817 cm<sup>-1</sup>, using a dilution analysis described by Filges et al. (Filges et al., 2018), who showed that the results obtained this way agree well with water vapor fractions measured with a conventional hygrometer. Figure 7 shows the measured amplitudes of the water and oxygen lines for samples of variable

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





humidity. Since the air came from a tank of constant composition, the oxygen concentration changes due to dilution of oxygen when water is added. Assuming that this is the sole cause of the change in measured absorption, since the line shapes were being constantly adjusted to account for changes in collisional broadening, it is straightforward to deduce the relation between the water fraction and the absorption amplitude. This calibration was used to generate the water fraction axes in Figures 4 and 6. We note that we did not take particular care to control or measure the quantity of dissolved gases, especially oxygen and carbon dioxide, in the water used for this experiment. While these gases would not significantly affect the water calibration, they may affect the water vapor correction of the oxygen measurement at the ppm level. More work needs to be done to investigate the water vapor correction of the oxygen measurement.

The observations described above were used to design a method to measure oxygen concentration in ambient air. Gas from the inlet to the analyzer is drawn through the cavity at a rate of about 100 scm and the conditions in the cavity are held stable at 340 hPa and 45° C. In its analysis mode the analyzer alternately measures ring-downs in the 7817 cm<sup>-1</sup> and 7878 cm<sup>-1</sup> regions. At 7878 cm<sup>-1</sup> measurements are made at 11 different frequencies, spaced by one FSR of the cavity and centered at the peak of the Q13Q13 line. Multiple ring-down measurements are made to improve the precision of the loss determination, with a total of 305 ring-downs allocated to one spectrum. In the 7817 cm<sup>-1</sup> region measurements are also made at 11 distinct frequencies at FSR spacings. Only 35 ring-downs are allocated to this spectral region, since the measurement of O<sub>2</sub> is much more important than water vapor. The data sets are analysed using a Levenberg-Marquardt fitting routine, which adjusts five free parameters in each region to find the best agreement to a spectral model based on Galatry line shapes, as described above. One of the outputs of the 7878 cm<sup>-1</sup> fit is the frequency offset of the FSR

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.



287

288

289

290

291

292

293

294

295

296

297

298

299

300

301

302

303

304

305

306

307

308

309

310



grid from the center of the Q13Q13 line. This information is used to adjust the position of the PZT actuated mirror to keep the measurements centered on the line, effectively stabilizing the optical path length of the cavity to the frequency of the O<sub>2</sub> line. The reported water fraction is obtained by multiplying the fitted amplitude of the water line by a calibration constant derived from the dilution experiment as explained above. For the O2 fraction a slightly more complicated procedure is followed. It was observed that the least-squares fitting of the data gives highly correlated results for the amplitude of the absorption line and the line width parameter v. The correlation may be due in part to the fitting procedure itself (Press et al., 1992) and it may also have a contribution from pressure variations that the pressure sensor is unable to detect. The ratio A<sub>O2</sub>/y can be determined from the fit much more precisely than A<sub>02</sub> alone and so gives a more sensitive measurement of molecular absorption. It also has the advantage of being independent of sample pressure, to the extent that the Galatry model applies (Figure 2). However, using the ratio A<sub>O2</sub>/y as a metric for absorption adds additional complications if measurements are to be made over a range of O<sub>2</sub> and water concentrations, because the O2/ N2 ratio and water concentration affect the line width independently of pressure and O2 concentration alone. To minimize systematic errors due to these broadening effects, we define a nominal y-parameter based on the measured amplitudes of the O2 and water lines and the line broadening dependences shown in Figures 3 and 4. The measured ratio A<sub>02</sub>/y is normalized by the nominal y to obtain a quantity that is ideally independent of pressure and water concentration, and this is the quantity that is multiplied by a calibration constant to give the reported O<sub>2</sub> fraction. In addition, a dry mole fraction is reported for O<sub>2</sub>, defined as the directly measured mole fraction corrected for water dilution.

The main goal in developing this instrument was to make high precision measurements of  $O_2$  mole fraction, based on absorption by the dominant  $^{16}O_2$  isotopologue.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.



311

312

313

314

315

316

317

318

319

320

321

322

323

324

325

326

327

328

329

330

331

333

334



The absorption lines of the rarer isotopologues are also present nearby, so a mode of operation was included in which one laser is scanned over neighboring lines of <sup>16</sup>O<sub>2</sub> and <sup>16</sup>O<sup>18</sup>O and the ratio of amplitudes is used to derive an isotopic ratio, reported in the usual delta notation. In this case the operating pressure was reduced to 160 hPa to improve the resolution of the nearby lines. The lines measured were the Q3Q3 line of <sup>16</sup>O<sub>2</sub>, at 7882.18670 cm<sup>-1</sup>, and the Q9Q9 line of <sup>16</sup>O<sup>18</sup>O, at 7882.050155 cm<sup>-1</sup>. The measurement procedure is very much like that for the O<sub>2</sub> fraction measurement, so it will not be described in detail, only the main differences will be noted. One is that in determining an isotopic ratio there is no advantage to be obtained from normalizing absorption amplitudes to line widths, instead we simply take the ratio of amplitudes to compute delta. Although the Q9Q9 line and its neighbor Q8Q8 are the strongest ones in this band, absorption by <sup>16</sup>O<sup>18</sup>O is still very weak, only about 5x10<sup>-9</sup> cm<sup>-1</sup> at the line center under the conditions we used. Consequently, the signal-to-noise that can be achieved with this analyzer is not adequate to determine both the amplitude and the width of the <sup>16</sup>O<sup>18</sup>O line with useful precision, so in the fitting step the v-parameter of the <sup>16</sup>O<sup>18</sup>O line is constrained to be a constant factor times the fitted y-parameter for the <sup>16</sup>O<sub>2</sub> line. Additionally, because of the weakness of the rare isotopologue absorption, the majority of ring-downs in each spectrum is devoted to measuring <sup>16</sup>O<sup>18</sup>O i.e. 232 ring-downs in each spectrum versus only 40 for <sup>16</sup>O<sub>2</sub>. This implies that the mole fraction measurement in the isotopic mode is much less precise than when the analyzer measures the Q13Q13 line alone.

### 3. Results and Discussions

#### 3.1. Laboratory tests at Picarro, Santa Clara

### 3.1.1. Temperature and pressure sensitivity

One set of tests was done to determine how well the goal was met of minimizing the susceptibility of the concentration measurements to uncontrolled noise or drift of the sample

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





temperature and pressure. For these tests the analyzer sampled dry synthetic air from a tank and the temperature and pressure setpoints of the cavity were adjusted upward and downward from the nominal values, to obtain an estimate of the differential response. We express the sensitivity to experimental conditions in relative form, that is the derivative with respect to temperature or pressure divided by the signal under nominal conditions.

From these experiments, we determined a temperature sensitivity of -2.1x10<sup>-4</sup> K<sup>-1</sup> and a pressure sensitivity of +9.8x10<sup>-6</sup> hPa<sup>-1</sup>. The temperature sensitivity is somewhat larger than expected based on a calculation using Hitran data to estimate the temperature dependences of all the quantities that go into the measured absorption of the Q13Q13 line. The pressure sensitivity is strikingly small, indicating a good cancelation of the pressure dependence of absorption amplitude and line width. Both temperature and pressure sensitivities are small enough to have a negligible effect on short-term precision of measurements made with the stabilized ring-down cavity, though long-term drifts in the sensors are always a matter of concern.

## 3.1.2. Measurement precision and Drift

Measurement precision was evaluated by analyzing synthetic air containing nominal atmospheric concentrations of CO<sub>2</sub> and CH<sub>4</sub> from an aluminum Luxfer cylinder over a period of several days. The tank, oriented horizontally and thermally insolated (though not controlled), was connected directly to the instrument (S/N TADS2001) with a 2-stage regulator and stainless-steel tubing and reducing the flow with an additional orifice to about 55 sccm. For the isotopic mode of operation, the precision of the measurement was also tested by making repeated measurements from a tank of clean, dry synthetic air.

Figure 8 shows the time series of the precision test data, displaying the reported oxygen concentration, the height of the oxygen absorption peak, the width of the oxygen

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





absorption peak and the ambient temperature. The residual drift of the analyzer, although small, is nevertheless significant given the stringent targets set forth by the WMO-GAW program. Possible sources of drift include: temperature drifts due to sensor drift or gradients; pressure errors due to sensor drift; optical artifacts such as parasitic reflections, higher order cavity mode excitation, and/or loss nonlinearity that can distort the reported oxygen spectrum. More work is required to identify and eliminate these small drifts.

The Allan standard deviation of the reported  $O_2$  fraction is shown in Figure 9. The ordinate on this plot is the square root of the Allan variance of reported mole fraction, so 1 ppm in these units corresponds to about 5 per meg in the ratio of  $O_2$ /  $N_2$ . The precision of averaged measurements improves as  $\tau^{-1/2}$  for approximately 5000 s and reaches 1 ppm in less than 10 minutes and remains below 1 ppm for time scales on the order of about 1 hour.

Figure 10 shows the precision of  $\delta(^{18}O)$  (uncalibrated) derived from the ratio of lines measured at 7882 cm<sup>-1</sup>. Because of the weak signal from the  $^{16}O^{18}O$  line, it is necessary to average for more than 20 seconds or more to achieve 1‰ precision on the isotopic ratio. As for the concentration measurement, averaging improves the measurement precision for times scale up to about 1 hour.

### 3.2. Laboratory measurements at the University of Bern

#### 3.2.1. Measurements of standard gases

The performance of the instrument was tested by analyzing eight standard gases with precisely known CO<sub>2</sub> and O<sub>2</sub> compositions (Table 1) using the CRDS analyzer and comparing it to parallel measurements with a paramagnetic oxygen sensor (PM1155 oxygen transducer, Servomex Ltd, UK) embedded to a commercially available Oxzilla fuel cell oxygen analyzer (OXZILLA II, Sable Systems International, USA) as well as with an isotope ratio mass spectrometer (IRMS, Finnigan Delta<sup>Plus</sup>XP). The design of the measurement set-up is shown

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.



383

384

385

386

387

388

389

390

391

392

393

394

395

396

397

398

399

400

401

402

403

404

405

406



multi-port valve (V2) was used to select among the standard gases. The flow from each cylinder was adjusted to about 120 ml min<sup>-1</sup> which was eventually directed to a selection valve (V1), allowing switching between ambient air and standard gases. Flow towards and out of the Oxzilla was controlled by the pressure controlling unit. The O2 mixing ratio of this incoming gas was first measured on the Paramagnetic O<sub>2</sub> sensor and then directed towards a non-dispersive infrared analyzer (NDIR) (Li-7000, LICOR, USA) for measuring CO2 and H<sub>2</sub>O. The outflow from this analyzer (100 ml min<sup>-1</sup>) returns to the pressure controlling unit and was eventually divided between the CRDS analyzer (which uses about 75-80 ml min<sup>-1</sup>) and the IRMS (~ 20 ml min<sup>-1</sup>) via a Tee-junction. Each cylinder was measured for two hours in each system controlled by a Lab VIEW program. In priori, we investigated the influence of this Tee-junction, which splits the gas flow between the CRDS and the IRMS, on the measured O2 values. Manning (2001) showed that the fractionation of O<sub>2</sub> in the presence of a Tee-Junction is strongly dependent on the splitting ratios as well as temperature and pressure gradients. Hence, we measured and compared the O<sub>2</sub> mixing ratios of two standard gases (CA07045 and CA060943) in two cases: i) in the presence of a Tee-junction with different CRDS to IRMS splitting ratios and ii) without a Tee-junction so that all gas flow is directed towards the CRDS analyzer. The splitting ratios in these test experiments vary from 1:1 to 1:100, and reversed to change the major flow direction either to the CRDS or the IRMS. Note that the experimental condition in this manuscript is with a 4:1 splitting ratio (i.e. ~ 80 ml min<sup>-1</sup> towards the CRDS analyzer and ~ 20 ml min<sup>-1</sup> towards the IRMS).

In the cases of the smaller splitting ratios (1:1, 1:4 and 4:1), which are relevant for the

results presented in this study, only minor differences in the measured O<sub>2</sub> mixing ratios were

in Figure 11. Standard gases were directly connected to the pressure controlling unit, and a

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





observed when compared to case b (i.e. without a Tee-junction). For these two cylinders measured, the average differences in these cases were about 0.5 ppm, calculated as the mean of the differences in the raw  $O_2$  measurements of the last 60 seconds. The negligible fractionation can indeed be the result of smaller splitting ratios while strong influence is usually expected in case of larger splitting ratios (Stephens et al., 2007). For higher splitting ratios, the result seems inconclusive without any dependence on the ratios due to the strong decline in the cylinder temperature (specifically at the pressure gauge) caused by higher flow to achieve the higher splitting ratios (as high as 1:100). Hence, these tests need to be conducted in a temperature controlled condition and the results could not be discussed in this manuscript.

Figure 12 shows the standard gas measurements for the seven cylinders with known CO<sub>2</sub> and O<sub>2</sub> mixing ratios (Table 1) using both the CRDS and the Paramagnetic analyzers. Standard eight, which has too high O<sub>2</sub>, is not shown in the figure as the figure is zoomed-in to better illustrate the change in O<sub>2</sub> for the remaining cylinders. While the first five cylinders contain O<sub>2</sub> and CO<sub>2</sub> fractions comparable to ambient air values, standards 6 & 8 had either very low and very high O<sub>2</sub>, respectively. In addition, standard 6 and 7 have very low and very high CO<sub>2</sub> mixing ratios. Note that due to its very high CO<sub>2</sub> content (~ 2700 ppm), standard 7 was not measured on the IRMS and hence the O<sub>2</sub> mixing ratios are unknown. The measured mixing ratios for the six standard gases between the two systems are in very good agreement while cylinder 7 showed an opposing signal for the two analyzers compared to standard 6 (Figure 12). While the Paramagnetic analyzer showed a higher O<sub>2</sub> mixing ratio, the values from the CRDS analyzer are lower in O<sub>2</sub>. This can be associated with the very high CO<sub>2</sub> mixing ratio in standard 7, which leads to a strong dilution effect in the CRDS analyzer as it does not include any correction function for dilution effect from CO<sub>2</sub>. However, such high

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





CO<sub>2</sub> mixing ratios may not be that important for most atmospheric research. Yet, it should be considered to include a parallel CO<sub>2</sub> mixing ratios measurement to the instrument as it will further improve the accuracy. This would be especially important for biological or physiological studies where a wide range of CO<sub>2</sub> and O<sub>2</sub> concentrations must be expected.

The measurement precision of the CRDS analyzer was calculated as the standard error of the mean i.e. the standard deviation  $(1-\sigma)$  of the last 1-minute raw measurements divided by the square root of the number of measurements (n=60), and for all these cylinders the values are usually between 0.5 ppm to 0.7 ppm. For parallel measurements of these cylinders using a Paramagnetic analyzer, we obtained a precision of about 1 ppm, calculated exactly the same way.

We also made a correlation plot to see which of the two instruments are in better agreement with the assigned values based on IRMS measurements for the individual cylinders. While similar correlation coefficients were observed for both analyzers, different slopes were calculated (Fig. A.1). This is due to the fact that the IRMS measures the  $O_2$  to  $N_2$  ratio ( $\delta(O_2/N_2)$ ) in per meg, while the CRDS and the Paramagnetic analyzers provide non-calibrated  $O_2$  mixing ratios in units of ppm and per meg, respectively. If we exclude the two standard gases with the highest and lowest  $O_2$  mixing ratios (standards 7 and 8) that are subjected to strong dilution effects, both the slope and the  $r^2$  values decrease from those shown in Figure A.1. But this decrease is larger in the case of the Paramagnetic measurements, implying a slightly better linearity of the CRDS analyzer.

#### 3.2.2. Measurements of ambient air

Ambient air measurements were conducted from the roof top of our laboratory at the University of Bern to evaluate the analyzer's performance under atmospheric variability. Ambient air was continuously aspirated from the inlet at the roof of the building at a flow rate

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.



455

456

457

458

459

460

461

462

463

464

465

466

467

468

469

470

471

472

473

474

475

476

477

478



of ~ 250 ml min<sup>-1</sup> which is then dried using a cooling trap kept at -90 °C towards the switching valve (V1) and measured in similar way to the standard gases as explained above. The measurement values obtained here were compared with the parallel measurements by the Paramagnetic sensor to test the instruments stability and accuracy.

Figures 13 panels a &b show the 1-minute average ambient air measurements from the rooftop inlet by the Paramagnetic and the CRDS analyzers at the beginning of the testing period including standard gases measured every 12-hour. While the Paramagnetic analyzer seems to be stable, the CRDS analyzer showed a strong drift for an extended period. This can be due to unstable conditions in the CRDS measurement system as it started operating right after it was unpacked. Hence, we looked into its DAS temperature and pressure records, which were stable within the manufacturer's recommended range during this period. As the CRDS analyzer incorporates a water correction function, interference from this species should be well accounted. Even comparing the analyzer's parallel water measurements to water measurements by the NDIR system such a drift was not observed. It should be noted that the two internal standard gases which were less frequently measured (every 12 hours) during this period were also drifting in similar pattern. This implies that the drift is associated with the analyzer. Interestingly, we observed that the two cylinders follow exactly the same drift pattern that can be modeled using a polynomial function which can then be used to correct for the observed drift in the ambient air measurements. After applying a polynomial drift correction, we were able to fully account for the observed drift. However, the manufacturer decided to further investigate possible causes of this drift. After further improvements, we obtained the first commercial analyzer in September 2017 and repeated the above tests (Figure 13 c &d). No such drift was observed any more in the standard gases or in ambient air measurements.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.



479

486

487

488

489

490

491

492

493

494

495

496

497

498

499

500

501



3.2.3. Water correction test

Measurements of oxygen are reported as both wet  $(O_{2, raw})$  and dry  $(O_{2, dry})$  mole fractions by the CRDS analyzer as it also measures water vapor in parallel at its water absorption line (7817 cm<sup>-1</sup>), and corrects for the dilution effect based on an inbuilt numerical function:

484 
$$O_{2,dry} = \frac{O_{2,raw}}{1 - f_{H,O}}$$
 (4)

where  $f_{H2O}$  is the measured water mole fraction.

The efficiency of water correction by this function was assessed in two ways: i) by comparing the water vapor content in standard air measured by this analyzer with similar measurements by the NDIR analyzer and ii) by comparing the oxygen mixing ratios between non-dried ambient air measured and corrected for water dilution by the CRDS analyzer with dried air measured using a paramagnetic analyzer.

Figure 14 shows the water vapor content for standard gases measured continuously for two days by the CRDS and the NDIR analyzers. Note that the two data sets are manually fitted to each other as the measured water values by the NDIR analyzer are not calibrated. Based on these plots, the two analyzers are in very good agreement although there are small differences during very dry conditions (low water content).

The water correction test was conducted by measuring dried ambient air (Figure 15a) into both analyzers as well as allowing non-dried air to the CRDS analyzer only (Figure 15b) and comparing the difference in O<sub>2</sub> measurements in both cases (Figures 15c & 15d). shows the water contents of dried ambient air measured in both analyzers while Figure 15b shows in case non-dried air is admitted to the picaro analyzer only (note that the CRDS uses its in-built water correction function). The measurements of the Paramagnetic analyzer were scaled to

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





ppm units by applying the correlation equation obtained from the six standard gas measurements of the two analyzers (Fig. A.1). Note that the CRDS measurements were corrected for the observed drift using the polynomial fit to the two standard gas measurements stated above.

In the first period of the measurement when both analyzers measured dried ambient air, the absolute differences between the 1-minute averages measured over two days by the two analyzers were mostly within 15 ppm and symmetrically distributed around zero. However, when wet air was admitted to the CRDS analyzer and the in-built water correction was applied, a stronger variability was observed in the calculated differences. This implies stronger short term variability in the CRDS analyzer measurement values (as nothing was changed for the Paramagnetic measurement system) when wet samples were analyzed. The more negative values in the differences can also be associated with overestimation of the O<sub>2</sub> mixing ratios by the CRDS originating from an overestimated water correction. However, detailed evaluation of the analyzer's water correction function is beyond the scope of this study.

### 3.3. Field Measurements

After a series of tests at University of Bern, we conducted multiple field measurements at the High Altitude Research Station Jungfraujoch and the Beromünster tall tower sites in Switzerland described below.

3.3.1. Tests at the High Altitude Research station Jungfraujoch

The High Alpine research station Jungfraujoch is located on the northern ridge of the Swiss Alps (46° 33′ N, 7° 59′ E) at an elevation of 3580 m a.s.l. It is one of the global atmospheric watch (GAW) stations well-equipped for measurements of numerous species and aerosols. The site is above the planetary boundary layer most of the time due to its high

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





elevation (Henne et al., 2010; Zellweger et al., 2003). However, thermally uplifted air from the surrounding valleys during hot summer days or polluted air from the heavily industrialized northern Italy may reach at this site (Zellweger et al., 2003). The Division of Climate and Environmental Physics at the University of Bern has been monitoring CO<sub>2</sub> and O<sub>2</sub> mixing ratios at this site based on weekly flask sampling and continuous measurements since 2000 and 2004, respectively. The CO<sub>2</sub> mixing ratio is measured using a commercial NDIR analyzer (S710 UNOR, SICK MAIHAK) while O<sub>2</sub> is measured using the Paramagnetic sensor (PM1155 oxygen transducer, Servomex Ltd, UK) and fuel cells (Max-250, Maxtec, USA) embedded within a home-built controlling unit. Similar to the comparison tests at the University of Bern, we have conducted parallel measurements between the CRDS analyzer and the paramagnetic cell at this high altitude site during 03 – 14 February 2017. The measurement of ambient air at the Jungfraujoch system is composed of sequential switching between a low span (LS) and high span (HS) calibration gases followed by a target gas (T) measurement (once a day) to evaluate the overall system performance and finally a working gas (WG) measurement before switching back to ambient air.

Figure 16 (top panel) shows the calibrated 1-minute averaged O<sub>2</sub> mixing ratios measured at this high altitude site in comparison with the Paramagnetic oxygen analyzer already available at the site. While a strong variability was observed during the measurement period of 10-days by both analyzers, a very good agreement was observed between them.

Figure 16 (bottom panel) shows the absolute difference of 1-minute averages in atmospheric  $O_2$  measured at Jungfraujoch between the two analyzers which are mostly within  $\pm 5$  ppm range (but sometimes going as high as  $\pm 10$  ppm) without an offset. However, for generally reported 10-minutes, half-hourly or hourly means these values correspond to < 1.5 ppm, < 1 ppm and < 0.65 ppm.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.



550

551

552

553

554

555

556

557

558

559

560

561

562

563

564

565

566

567

568

569

570

571

572

573



3.3.2. Tests at the Beromünster tall tower site

The Beromünster tower is located near the southern border of the Swiss Plateau, the comparatively flat part of Switzerland between the Alps in the south and the Jura mountains in the northwest (47° 11′ 23" N, 8° 10' 32" E, 797 m a.s.l.), which is characterized by intense agriculture and rather high population density. A detailed description of the tower measurement system as well as a characterization of the site with respect to local meteorological conditions, seasonal and diurnal variations of greenhouse gases, and regional representativeness can be obtained from previous publications (Berhanu et al., 2016; Berhanu et al., 2017; Oney et al., 2015; Satar et al., 2016). The tower is 217.5 m tall with access to five sampling heights (12.5 m, 44.6 m, 71.5 m, 131.6 m, 212.5 m) for measuring CO, CO<sub>2</sub>, CH<sub>4</sub> and H<sub>2</sub>O using Cavity Ring Down Spectroscopy (Picarro Inc., G-2401). By sequentially switching from the highest to the lowest level, mixing ratios of these trace gases were recorded continuously for three minutes per height, but only the last 60 seconds were retained for data analysis. The calibration procedure for ambient air includes measurements of reference gases with high and low mixing ratios traceable to international standards (WMO-X2007 for CO<sub>2</sub> and WMO-X2004 for CO and CH<sub>4</sub>), as well as target gas and more frequent working gas determinations to ensure the quality of the measurement system. From two years of data a long-term reproducibility of 2.79 ppb, 0.05 ppm, and 0.29 ppb for CO, CO<sub>2</sub> and CH<sub>4</sub>, respectively was determined for this system (Berhanu et al., 2016).

Between 15.02.2017 and 02.03.2017, we have connected the new CRDS oxygen analyzer in series with the CO<sub>2</sub> analyzer (Picarro G-2401) and measured the O<sub>2</sub> mixing ratios at the corresponding heights. Similar to the CO<sub>2</sub> measurements, O<sub>2</sub> was also measured for three minutes at each height. During this period, we have evaluated the two features (isotopic mode and concentration mode) of the CRDS analyzer. In the isotopic mode, the CRDS

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





measures the  $\delta^{18}O$  values as well as the  $O_2$  concentration while in concentration mode only the latter was measured.

During the tests conducted at this tower site, we first evaluated the two operational modes (concentration vs isotopic modes) of the CRDS analyzer. Ambient air measurements on isotopic mode over a 4-days period showed a strong variability in the measured oxygen mixing ratios and it was not possible to distinguish the variability in the  $O_2$  mixing ratios among the five height levels. The calculated 1-minute standard error for ambient air measurements was as high as 10 ppm while a standard error of less than 1 ppm was determined from similar measurements in the concentration mode. Additionally, comparing the  $O_2$  values between the two modes, frequent short time variation in ambient air  $O_2$  (~ 200 ppm) was observed in the isotope mode measurements while the variation in the concentration mode is significantly smaller (~ 30 ppm). This precision degradation is due to the weaker  $^{16}O$  oxygen line used for the isotopic mode, and the fact that far more ring-downs are collected on the rare isotopologue in isotopic mode Hence, we have conducted the remaining test measurements in concentration mode.

As this tower has five sampling height levels, we first followed three minutes of switching per inlet level, which enables four measurements per hour at a given level. However, we noticed hardly any difference among the different levels due to strong short term variability in O<sub>2</sub> mixing ratios between the consecutive heights. Hence, we switched to a longer sampling period of six-minutes per height. Figure 17 shows the diurnal CO<sub>2</sub> and O<sub>2</sub> variations at the lowest (12 m) and highest (212.5 m) sampling heights of the tower. These two heights were selected simply to better illustrate the difference in the mixing ratios. The CO<sub>2</sub> mixing ratios on the top panel show higher values at the 12 m inlet than the highest level most of the day due to its closeness to sources except during the afternoon (11:00 - 17:00

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





UTC) when both levels show similar but decreasing  $CO_2$  mixing ratios. This is due to presence of a well-mixed planetary boundary layer (PBL) (Satar et al., 2016). The lag in  $CO_2$  peak between the two height levels by about two hours indicates the duration for uniform vertical mixing along the tower during winter 2017. The opposite variability patterns are also clearly visible in the  $O_2$  mixing ratios shown in the lower panel with a clear distinction between the two height levels during early in the morning and in the evening while similar  $O_2$  values were observed in the afternoon. These opposing profiles are expected as  $CO_2$  and  $O_2$  are linearly coupled with a mean oxidation ratio of -1.1  $\pm$  0.05 (Severinghaus, 1995) for land-biospheric processes (photosynthesis and respiration) and -1.44  $\pm$  0.03 for fossil fuel burning (Keeling, 1988b).

Table 2 shows the oxidation ratios derived as the slopes of the linear regression between  $CO_2$  and  $O_2$  mixing ratios at the different height levels measured on 25 February 2017. Accordingly, height dependent slopes were observed with a slope of  $-0.98 \pm 0.06$  at the lowest level, close to the biological processes induced slope but slightly lower than its mean value. For the highest level, we calculated a slope of  $-1.60 \pm 0.07$  a value close to fossil fuel combustion oxidation ratio. Note that depending on fossil fuel type the oxidation ratio can range between -1.17 and -1.95 for coal and natural gas, respectively (Keeling, 1988b). While the slopes derived for the two other levels (44.6 m and 131.6 m) show similar values between the highest and lowest height levels, possibly from mixed sources, the middle level showed a slightly higher slope than these two levels but still in the large range between the lowest and highest inlet heights.

# 3.4. Evaluation of the $\delta^{18}$ O measurements

To further evaluate the analyzer's performance in measuring stable oxygen isotopes, we conducted ambient air isotopic composition measurements as well as analyzed a standard

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.



622

623

624

625

626

627

628

629

630

631

632

633

634

635

636

637

638

639

640

641

642

643

644

645



gas without  $CO_2$  which has a known  $\delta^{18}O$  value. The choice of this  $CO_2$ -free air standard gas is twofold: one it has a known  $\delta^{18}O$  value and second as it has no  $CO_2$  possible interference from band overlap is avoided. For this test three 0.5 L glass flasks were preconditioned and filled with this standard gas to ambient pressure. These flasks were attached before or after the water trap (Fig. 11) and measured similar to ambient air measurements. These measurements were then compared with  $\delta(^{34}O/^{32}O)$  values obtained by parallel measurements using our IRMS.

Figure 18 shows the  $\delta^{18}$ O values of ambient air from the roof top with three consecutive measurements of glass flasks filled with CO<sub>2</sub>-free air in-between followed by a fourth flask filled with breath air. An excellent agreement was observed for measurements from both instruments for the three flasks filled with a standard gas. However, the fourth flask with breath air showed a signal opposite to the measurements by the IRMS. As breath air contains large amount of water CO2 in addition to O2, which can possibly interfere with the CRDS analyzer measurements, we have removed H<sub>2</sub>O and CO<sub>2</sub> by using a cryogenic trap (-130 °C) and in an additional experiment using Schütze reagent to remove both CO and CO<sub>2</sub>. However, we have not observed any improvement towards an agreement with the IRMS measurements. Therefore, any other gas component in the breath air must be relevant for the interference. Based on the absorption lines in the spectral range of the instrument (7878 cm<sup>-1</sup>) retrieved from HITRAN database, we expect interference either from carbon monoxide (now excluded by the tests) or methane or VOCs including acetone, ethanol, methanol or isoprenes, all of which have been measured in breath air (Gao et al., 2017; Gottlieb et al., 2017; Mckay et al., 1985; Ryter and Choi, 2013; Wolf et al., 2017). Further investigations have to shed light on these interferences in order to take corresponding action to surpass these shortcomings in the isotope analysis based on cavity ring-down spectroscopy.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





#### 4. Conclusions

We have thoroughly evaluated the performance of a new CRDS analyzer which measures O<sub>2</sub> mixing ratios and isotopic composition combining laboratory and field tests. Even if a drift in the analyzer was observed at the beginning of this study, which if it appears can be easily corrected by calibration, the recent analyzers built by the manufacturer did not show such instrumental drift. However, prior tests are recommended to see the analyzer's stability.

The T-split tests for the current measurement setup based on the measurements of two standard gases showed a difference within the measurement uncertainty. However, this effect may become significant while applying larger splitting ratios and we recommend conducting further experiments to accurately quantify this influence for larger splitting ratios.

We have observed a strong influence of dilution in the measured O<sub>2</sub> values during the presence of high CO<sub>2</sub> mixing ratios. Even if such an influence may not be critical for the present study, such an effect might be significant in other studies where higher CO<sub>2</sub> mixing ratios might be present and we recommend following a correction strategy based on parallel CO<sub>2</sub> measurements. This also applies for more accurate analysis.

The water correction applied by the instrument's in-built function seems to sufficiently correct for the water vapor influence. However, a larger variability of the difference was observed between the CRDS analyzer and the Paramagnetic cell when dried samples were used in both systems. This can possibly be due to an overcorrection by the water correction function of the CRDS analyzer when dried samples were used. This is particularly true for the very low water vapor range (< 100 ppm).

Based on the analysis of  $O_2$  mixing ratios in the concentration and isotopic modes, we have observed about a significant decrease in precision (about ten-fold) in the latter

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





measurement mode. The measured  $\delta^{18}$ O values for the standard air by the CRDS analyzer are in excellent agreement with the IRMS values. However, such measurements for a breath air showed a contrasting signal, possibly due to interference from other gases as breath air contains CO<sub>2</sub>. CH<sub>4</sub> and CO in addition to oxygen. Hence, we recommend further investigation on such possible contaminants and how to possibly remove them while conducting ambient air measurements. Acknowledgement We would like to thank ICOS-RI and the Swiss National Science Foundation (SNF) for funding ICOS-CH (20FI21\_148994, 20FI21\_148992). We are also grateful to the International Foundation High Alpine Research Stations Jungfraujoch and Gornergrat. The measurement system at the Beromünster tower was built and maintained by the CarboCount-CH (CRSII2\_136273) and IsoCEP (200020\_172550) projects both funded by SNF.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





694 List of Tables

Table 1. Assigned mixing ratios of standard gases used in this study and their corresponding

values measured by the NDIR, CRDS and IRMS at the University of Bern. <sup>1</sup>The assigned

values are based on measurements from different institutions (University of Bern (UB),

698 Scripps or NOAA, see column cylinder name). <sup>2</sup>Measurements are on the Bern scale for CO<sub>2</sub>

and O<sub>2</sub>. The Bern scale is shifted by +550 per meg. <sup>3</sup>Values on the Scripps scale.

700

Cylinder name	Assigned CO <sub>2</sub> (ppm) <sup>1</sup>	Assigned O <sub>2</sub> (per meg) <sup>1</sup>	CO <sub>2</sub> -IRMS (ppm) <sup>2</sup>	CO <sub>2</sub> - NDIR (ppm) <sup>2</sup>	O <sub>2</sub> -IRMS (per meg) <sup>2</sup>	O <sub>2</sub> - Paramagnetic (per meg) <sup>2</sup>	O <sub>2</sub> - CRDS (per meg) <sup>2</sup>
ST-1 LUX3576- UB	427.47	-1026	427.47	427.59	-1026	-1070	-1057
ST-2 LK922131- UB	368.09	599	368.09	367.82	599	560	590
ST-3 CA07045- Scripps	382.303	-271.6	382.50	381.99	278 (-272.2) <sup>3</sup>	302	281
ST-4 CA07043- Scripps	390.528	-476.4	390.69	390.15	71 (-479.5) <sup>3</sup>	66	63
ST-5 CA07047- Scripps	374.480	-807.7	374.70	374.17	-253 (-803.3) <sup>3</sup>	-212	-233
ST-6 CA04556- NOAA	192.44	-3410	191.21	191.64	-3410	-2905	-3013
ST-7 CA06943-	2699.45	-		2612.80	-	-2691	-3369

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





NOAA							
ST-8	411.49	37794	411.49	406.25	37794	34513	36017
LK76852-							
UB							

703 Table 2. The CO<sub>2</sub> and O<sub>2</sub> correlation coefficients at the different height levels derived using

704 the least square fit and the correlation coefficients  $(r^2)$ . Uncertainties are calculated as

standard error of the slope.

Height	Oxidation Ratios
	(O <sub>2</sub> :CO <sub>2</sub> )
	(32.332)
10.5	0.00 . 0.06 (0.40)
12.5 m	$-0.98 \pm 0.06  (0.48)$
44.6 m	$-1.29 \pm 0.07 (0.50)$
	1.25 = 0.07 (0.00)
71.5 m	$-1.49 \pm 0.08 (0.47)$
131.6 m	$-1.23 \pm 0.05 (0.55)$
	` ′
212.5 m	$-1.60 \pm 0.07 (0.61)$
	` ′

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





# 714 List of Figures




Figure 1. The Q13Q13 line of  $O_2$  measured in a sample of synthetic air at a sample temperature and pressure of 45° C and 333 hPa, respectively.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.









Figure 2. Best-fit values for the Galatry parameters of the Q13Q13 line of  $O_2$ , as a function of pressure. The line broadening parameter y is represented by circles and the line narrowing parameter z by squares. The solid lines are linear fits to the measurements. The best-fit offset and slope are 0.0227 and 0.004082 hPa<sup>-1</sup> for y, and -0.0169 and 0.001424 hPa<sup>-1</sup> for z.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.







734 Figure 3. Galatry parameters of the Q13Q13 line of  $O_2$  at 340 hPa and 45° C as a function of

735  $O_2$  mole fraction in binary  $O_2$  -  $N_2$  mixtures.

733

736 The linear fits to the data are  $y = 1.417 - 0.023 \text{ x f}_{02}$  and  $z = 0.481 + 0.085 \text{ x f}_{02}$ .

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.



737

738

739






Figure 4. Galatry parameters of the  $7816.75210~cm^{-1}$  water line in air at 340~hPa and  $45^{\circ}$  C as a function of water mole fraction. Black points are from measurements and red lines are linear fits:  $y = 1.7846 + 8.01~x~f_{H2O}$  and  $z = 0.656 + 3.60~x~f_{H2O}$ .

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.








Figure 5. Upper panel: spectrum of water in nitrogen (points) and fit to Voigt model (blue curve). Lower panel: Oxygen (green) and water (blue) lines in the Hitran database.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.








Figure 6. Galatry collisional broadening parameter of the oxygen Q13Q13 line at 340 hPa and 45° C versus water mole fraction. Black points are from measurements and the red line is a linear fit:  $y = 1.4109 + 0.467 f_{H2O}$ .

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.









Figure 7. Measured absorption line amplitudes for oxygen and water vapor for water vapor mixing ratios ranging from nearly 0 to 0.025. Black points are from measurements and the red line is a linear fit: with intercept  $7.78001 \times 10^{-6} \text{ cm}^{-1}$  and slope -0.014807.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





762763

764

765

766


Figure 8. Time series from a measurement of a single tank over about a week. The four panels shown the water-corrected oxygen concentration, the absorption peak loss minus the baseline loss, the measured Lorentzian broadening factor, and the ambient temperature (measured in the instrument housing), respectively. A windowed average of 300 seconds was applied to all four data sets.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.







768769 Figure

Figure 9. Precision of  $O_2$  mole fraction measured from a tank of synthetic air. Filled circles are measurements and the line shows the ideal  $\tau^{-1/2}$  dependence.

771

770

772

773

774

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.









Figure 10. Precision of  $\delta(^{18}O)$  measured from a tank of synthetic air. Filled circles are measurements and the line shows the ideal  $\tau^{-1/2}$  dependence.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.







788 Figure 11. Schematics of the measurement system used to compare the Picarro analyzer with

789 the Mass Spectrometer at Bern.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.








Figure 12. Comparison of oxygen mixing ratios for the seven standard gases measured using the CRDS analyzer (black) and the Paramagnetic sensors (red).

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.








Figure 13. Parallel ambient air measurements by the Paramagnetic and CRDS analyzers at the beginning of the testing period (Panels a & b, January 2017) and the second phase of testing (Panels c & d, September 2017). The spikes are measurements from the two standard gases bracketing the ambient air values.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.








Figure 14. Parallel water vapor measurements for a dried ambient air by both the NDIR and CRDS analyzers. Note that the water values from the NDIR analyzer are not calibrated.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.








Figure 15. Results of water correction tests. Water measurements of the NDIR (left scale) for dry conditions (a,b) and the CRDS analyzer (right scale) for dry (a) and wet (b) conditions. The difference in oxygen measurements between the Paramagnetic and the CRDS instrument using the built-in water correction for the CRDS values under dry (c) and wet (d) conditions. Panels (e) and (f) show the population density functions.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.






Figure 16. Calibrated ambient air oxygen measurements (1-minute average) at the Jungfraujoch site using the CRDS and Paramagnetic analyzers both in ppm units (a) and the absolute difference between the two measurements in ppm (b) by matching time stamps.

868 869

864

865

866

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.








Figure 17. Diurnal variations of  $CO_2$  (top) and  $O_2$  (bottom) measurements from the 12 m (red) and the 212.5 m (black) height levels at Beromünster tower.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.








Figure 18. Consecutive  $\delta^{18}O$  measurements of a standard gas (CO<sub>2</sub>-free air) filled into three flasks followed by measurement of breath air using the CRDS analyzer (top) and IRMS (bottom). These measurements were carried out in the middle of ambient air measurements.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





## 896 Appendix A.

## Additional plots




Figure A.1. Correlations between the  $O_2$  mixing ratios measured by the CRDS and Paramagnetic analyzers with the mass spectrometric measurements (uncalibrated values). The left panels are for all the cylinders measured (standards 1 to 8) while the right ones are after zooming only to standards 1-5.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





## 911 References

- 912 Battle, M., Bender, M. L., Tans, P. P., White, J. W. C., Ellis, J. T., Conway, T., and Francey, R. J.: Global
- 913 carbon sinks and their variability inferred from atmospheric O-2 and delta C-13, Science, 287, 2467-
- 914 2470, 2000.
- 915 Bender, M. L., Tans, P. P., Ellis, J. T., Orchardo, J., and Habfast, K.: A High-Precision Isotope Ratio
- 916 Mass-Spectrometry Method for Measuring the O-2 N-2 Ratio of Air, Geochim Cosmochim Ac, 58,
- 917 4751-4758, 1994.
- 918 Berhanu, T. A., Satar, E., Schanda, R., Nyfeler, P., Moret, H., Brunner, D., Oney, B., and Leuenberger,
- 919 M.: Measurements of greenhouse gases at Beromünster tall tower station in Switzerland, Atmos.
- 920 Meas. Tech., 9, 2016.
- 921 Berhanu, T. A., Szidat, S., Brunner, D., Satar, E., Schanda, R., Nyfeler, P., Battaglia, M., Steinbacher,
- 922 M., Hammer, S., and Leuenberger, M.: Estimation of the fossil-fuel component in atmospheric CO2
- 923 based on radiocarbon measurements at the Beromünster tall tower, Switzerland, Atmos. Chem.
- 924 Phys. Discuss., 2017, 1-33, 2017.
- 925 Crosson, E. R. J. A. P. B.: A cavity ring-down analyzer for measuring atmospheric levels of methane,
- 926 carbon dioxide, and water vapor, 92, 403-408, 2008.
- 927 Filges, A., Gerbig, C., Rella, C. W., Hoffnagle, J., Smit, H., Krämer, M., Spelten, N., Rolf, C., Bozóki, Z.,
- 928 Buchholz, B., and Ebert, V.: Evaluation of the IAGOS-Core GHG Package H2O measurements during
- 929 the DENCHAR airborne inter-comparison campaign in 2011, Atmos. Meas. Tech. Discuss., doi:
- 930 10.5194/amt-2018-36, 2018. 2018.
- 931 Fleisher, A. J., Hodges, J., and Sironneau, V.: Collision-dependent line areas in the a1 $\Delta g \leftarrow X3\Sigma g$  band
- 932 of molecular oxygen, 2015.
- 933 Gao, F., Zhang, X., Zhang, X., Wang, M., and Wang, P.: Virtual electronic nose with diagnosis model
- 934 for the detection of hydrogen and methane in breath from gastrointestinal bacteria, 28-31 May 2017
- 935 2017, 1-3.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.



936



937 A., Chance, K., Drouin, J., Flaud, J., Gamache, R. R., Hodges, J., Jacquemart, D., Perevalov, I., Perrin, A., 938 Shine, P., Smith, M., Tennyson, J., Toon, G., Tran, H., Tyuterev, G., Barbe, A., Császár, G., Devi, M., 939 Furtenbacher, T., Harrison, J., Hartmann, J., Jolly, A., Johnson, J., Karman, T., Kleiner, I., Kyuberis, A. 940 A., Loos, J., Lyulin, M., Massie, S., Mikhailenko, S., Moazzen-Ahmadi, N., Muller, S., Naumenko, O. V., 941 Nikitin, A. V., Polyansky, O. L., Rey, M., Rotger, M., Sharpe, S., Sung, K., Starikova, E., Tashkun, S., 942 Auwera, J., Wagner, G., Wilzewski, J., Wcisło, P., Yu, S., and Zak, E. J.: The HITRAN2016 molecular 943 spectroscopic database, 203, 3 - 69, 2017. 944 Goto, D., Morimoto, S., Ishidoya, S., Aoki, S., and Nakazawa, T.: Terrestrial biospheric and oceanic 945 CO2 uptake estimated from long-term measurements of atmospheric CO2 mole fraction,  $\delta 13C$  and δ(O2/N2) at Ny-Ålesund, Svalbard, Journal of Geophysical Research: Biogeosciences, doi: 946 947 10.1002/2017JG003845, 2017. n/a-n/a, 2017. 948 Gottlieb, K., Le, C. X., Wacher, V., Sliman, J., Cruz, C., Porter, T., and Carter, S.: Selection of a cut-off 949 for high- and low-methane producers using a spot-methane breath test: results from a large north

Gordon, E., Rothman, S., Hill, C., Kochanov, V., Tan, Y., Bernath, P., Birk, M., Boudon, V., Campargue,

- 951 Rep, 5, 193-199, 2017.
- 952 Hartmann, J.-M., Boulet, C., and Robert, D.: Collisional Effects on Molecular Spectra, Elsevier Science,

American dataset of hydrogen, methane and carbon dioxide measurements in breath, Gastroenterol

953 2008.

- 954 Henne, S., Brunner, D., Folini, D., Solberg, S., Klausen, J., and Buchmann, B.: Assessment of
- 955 parameters describing representativeness of air quality in-situ measurement sites, Atmos. Chem.
- 956 Phys., 10, 3561-3581, 2010.
- 957 Hodges, J. T., Layer, H. P., Miller, W. W., and Scace, G. E.: Frequency-stabilized single-mode cavity
- 958 ring-down apparatus for high-resolution absorption spectroscopy, 75, 849-863, 2004.
- 959 Keeling, R. F.: Development of an Interferometric Oxygen Analyzer for Precise Measurement of the
- 960 Atmospheric O2 Mole Fraction, UMI, 1988a.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





- 961 Keeling, R. F.: Measuring correlations between atmospheric oxygen and carbon dioxide mole
- 962 fractions: A preliminary study in urban air, J Atmos Chem, 7, 153-176, 1988b.
- 963 Keeling, R. F. and Manning, A. C.: 5.15 Studies of Recent Changes in Atmospheric O2 Content A2 -
- 964 Holland, Heinrich D. In: Treatise on Geochemistry (Second Edition), Turekian, K. K. (Ed.), Elsevier,
- 965 Oxford, 2014.
- 966 Keeling, R. F. and Shertz, S. R.: Seasonal and Interannual Variations in Atmospheric Oxygen and
- 967 Implications for the Global Carbon-Cycle, Nature, 358, 723-727, 1992.
- 968 Keeling, R. F., Stephens, B. B., Najjar, R. G., Doney, S. C., Archer, D., and Heimann, M.: Seasonal
- 969 variations in the atmospheric O-2/N-2 ratio in relation to the kinetics of air-sea gas exchange, Global
- 970 Biogeochem Cy, 12, 141-163, 1998.
- 971 Lamouroux, J., Sironneau, V., Hodges, J. T., and Hartmann, J. M.: Isolated line shapes of molecular
- 972 oxygen: Requantized classical molecular dynamics calculations versus measurements, Physical
- 973 Review A, 89, 042504, 2014.
- 974 Le Quéré, C., Andrew, R. M., Friedlingstein, P., Sitch, S., Pongratz, J., Manning, A. C., Korsbakken, J. I.,
- 975 Peters, G. P., Canadell, J. G., Jackson, R. B., Boden, T. A., Tans, P. P., Andrews, O. D., Arora, V. K.,
- 976 Bakker, D. C. E., Barbero, L., Becker, M., Betts, R. A., Bopp, L., Chevallier, F., Chini, L. P., Ciais, P.,
- 977 Cosca, C. E., Cross, J., Currie, K., Gasser, T., Harris, I., Hauck, J., Haverd, V., Houghton, R. A., Hunt, C.
- 978 W., Hurtt, G., Ilyina, T., Jain, A. K., Kato, E., Kautz, M., Keeling, R. F., Klein Goldewijk, K., Körtzinger, A.,
- 979 Landschützer, P., Lefèvre, N., Lenton, A., Lienert, S., Lima, I., Lombardozzi, D., Metzl, N., Millero, F.,
- 980 Monteiro, P. M. S., Munro, D. R., Nabel, J. E. M. S., Nakaoka, S. I., Nojiri, Y., Padín, X. A., Peregon, A.,
- 981 Pfeil, B., Pierrot, D., Poulter, B., Rehder, G., Reimer, J., Rödenbeck, C., Schwinger, J., Séférian, R.,
- 982 Skjelvan, I., Stocker, B. D., Tian, H., Tilbrook, B., van der Laan-Luijkx, I. T., van der Werf, G. R., van
- 983 Heuven, S., Viovy, N., Vuichard, N., Walker, A. P., Watson, A. J., Wiltshire, A. J., Zaehle, S., and Zhu,
- D.: Global Carbon Budget 2017, Earth Syst. Sci. Data Discuss., 2017, 1-79, 2017.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.



1008



985 Manning, A.: Temporal variability of atmospheric oxygen from both continuos and measurements 986 and a flask sampling network: tools for studying the global carbon cycle, Ph.D. Ph.D., University of 987 California, San Diego, San Diego, Calofornia, USA, 2001. 988 Manning, A. C. and Keeling, R. F.: Global oceanic and land biotic carbon sinks from the Scripps 989 atmospheric oxygen flask sampling network, Tellus B, 58, 95-116, 2006. 990 Manning, A. C., Keeling, R. F., and Severinghaus, J. P.: Precise atmospheric oxygen measurements 991 with a paramagnetic oxygen analyzer, Global Biogeochem Cy, 13, 1107-1115, 1999. 992 Marrero, T. R. and Mason, E. A.: Gaseous Diffusion Coefficients, Journal of Physical and Chemical 993 Reference Data 1, 3, 1972. 994 Martin, N. A., Ferracci, V., Cassidy, N., and Hoffnagle, J. A. J. A. P. B.: The application of a cavity ring-995 down spectrometer to measurements of ambient ammonia using traceable primary standard gas mixtures, 122, 219, 2016. 996 997 Mckay, L. F., Eastwood, M. A., and Brydon, W. G.: Methane Excretion in Man - a Study of Breath, 998 Flatus, and Feces, Gut, 26, 69-74, 1985. 999 Nevison, C. D., Keeling, R. F., Kahru, M., Manizza, M., Mitchell, B. G., and Cassar, N.: Estimating net 1000 community production in the Southern Ocean based on atmospheric potential oxygen and satellite 1001 ocean color data, Global Biogeochem Cy, 26, 2012. 1002 Oney, B., Henne, S., Gruber, N., Leuenberger, M., Bamberger, I., Eugster, W., and Brunner, D.: The 1003 CarboCount CH sites: characterization of a dense greenhouse gas observation network, Atmos. 1004 Chem. Phys., 15, 11147-11164, 2015. 1005 Press, W. H., Teukolsky, S. A., Vetterling, W. T., and Flannery, B. P.: Numerical Recipes 3rd Edition: 1006 The Art of Scientific Computing, Cambridge Printing Press, Cambridge, England, 1986. 1007 Press, W. H., Teukolsky, S. A., Vetterling, W. T., and Flannery, B. P.: Numerical recipes in C: the art of

scientific computing, Cambridge University Press, London, 1992.

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.



1032



1009 Rothman, L. S., Gordon, I. E., Babikov, Y., Barbe, A., Chris Benner, D., Bernath, P. F., Birk, M., 1010 Bizzocchi, L., Boudon, V., Brown, L. R., Campargue, A., Chance, K., Cohen, E. A., Coudert, L. H., Devi, V. 1011 M., Drouin, B. J., Fayt, A., Flaud, J. M., Gamache, R. R., Harrison, J. J., Hartmann, J. M., Hill, C., Hodges, 1012 J. T., Jacquemart, D., Jolly, A., Lamouroux, J., Le Roy, R. J., Li, G., Long, D. A., Lyulin, O. M., Mackie, C. 1013 J., Massie, S. T., Mikhailenko, S., Müller, H. S. P., Naumenko, O. V., Nikitin, A. V., Orphal, J., Perevalov, V., Perrin, A., Polovtseva, E. R., Richard, C., Smith, M. A. H., Starikova, E., Sung, K., Tashkun, S., 1014 1015 Tennyson, J., Toon, G. C., Tyuterev, V. G., and Wagner, G.: The HITRAN2012 molecular spectroscopic 1016 database, Journal of Quantitative Spectroscopy and Radiative Transfer, 130, 4-50, 2013. 1017 Ryter, S. W. and Choi, A. M. K.: Carbon monoxide in exhaled breath testing and therapeutics, J Breath 1018 Res, 7, 2013. 1019 Satar, E., Berhanu, T. A., Brunner, D., Henne, S., and Leuenberger, M.: Continuous CO2/CH4/CO 1020 measurements (2012-2014) at Beromünster tall tower station in Switzerland, Biogeosciences, 13, 1021 2623-2635, 2016. 1022 Severinghaus, J. P.: Studies of the terresterial O2 and carbon cycles in sand dune gases and in 1023 Biosphere Doctoral Ph.D., Columbia University, New York, USA, 1995. 1024 Steig, E. J., Gkinis, V., Schauer, A. J., Schoenemann, S. W., Samek, K., Hoffnagle, J., Dennis, K. J., and 1025 Tan, S. M.: Calibrated high-precision <sup>17</sup>O-excess measurements using cavity ring-down 1026 spectroscopy with laser-current-tuned cavity resonance, Atmos. Meas. Tech., 7, 2014. 1027 Stephens, B. B., Bakwin, P. S., Tans, P. P., Teclaw, R. M., and Baumann, D. D.: Application of a 1028 differential fuel-cell analyzer for measuring atmospheric oxygen variations, J Atmos Ocean Tech, 24, 1029 82-94, 2007. 1030 Tennyson, J., Bernath, P. F., Campargue, A., Császár, A. G., Daumont, L., Gamache, R. R., Hodges, J. T., 1031 Lisak, D., Naumenko, O. V., Rothman, L. S., Tran, H., Zobov, N. F., Buldyreva, J., Boone, C. D., De Vizia,

M. D., Gianfrani, L., Hartmann, J.-M., McPheat, R., Weidmann, D., Murray, J., Ngo, N. H., and

Manuscript under review for journal Atmos. Meas. Tech.

Discussion started: 17 January 2019 © Author(s) 2019. CC BY 4.0 License.





1033	Polyansky, O. L.: Recommended isolated-line profile for representing high-resolution spectroscopic
1034	transitions (IUPAC Technical Report), 86, 1931–1943, 2014.
1035	Tohjima, Y.: Method for measuring changes in the atmospheric O-2/N-2 ratio by a gas
1036	chromatograph equipped with a thermal conductivity detector, J Geophys Res-Atmos, 105, 14575-
1037	14584, 2000.
1038	Valentino, F. L., Leuenberger, M., Uglietti, C., and Sturm, P.: Measurements and trend analysis of O2,
1039	CO2 and $\delta$ 13C of CO2 from the high altitude research station Junfgraujoch, Switzerland — A
1040	comparison with the observations from the remote site Puy de Dôme, France, Sci Total Environ, 391,
1041	203-210, 2008.
1042	Varghese, P. L. and Hanson, R. K.: Collisional narrowing effects on spectral line shapes measured at
1043	high resolution, Appl. Opt., 23, 2376-2385, 1984.
1044	Wójtewicz, S., Cygan, A., Masłowski, P., Domysławska, J., Wcisło, P., Zaborowski, M., Lisak, D.,
1045	Trawiński, R. S., and Ciuryło, R.: Spectral line-shapes of oxygen B-band transitions measured with
1046	cavity ring-down spectroscopy, Journal of Physics: Conference Series, 548, 012028, 2014.
1047	Wolf, P. G., Parthasarathy, G., Chen, J., O'Connor, H. M., Chia, N., Bharucha, A. E., and Gaskins, H. R.:
1048	Assessing the colonic microbiome, hydrogenogenic and hydrogenotrophic genes, transit and breath
1049	methane in constipation, Neurogastroent Motil, 29, 2017.
1050	Zellweger, C., Forrer, J., Hofer, P., Nyeki, S., Schwarzenbach, B., Weingartner, E., Ammann, M., and
1051	Baltensperger, U.: Partitioning of reactive nitrogen (NO <sub>y</sub> ) and dependence on
1052	meteorological conditions in the lower free troposphere, Atmos. Chem. Phys., 3, 779-796, 2003.